DECLARACIÓN DE CONFORMIDAD CE
CE DECLARATION OF CONFORMITY
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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- ◆ Lenses must not be disassembled by the user. Were there any dust or dirt to be cleaned, clean it with a natural horse hair brush or a smooth piece of cloth, fluff-free, dampened with a bit of xylene or toluene.
- ◆ To remove dust from lenses, blow them with a plastic bulb or clean them with a natural horsehair brush.
- ◆ Use non-corrosive lubricants in metallic parts, being careful of not touching optical parts.

NOTE: Never disassemble the objective lenses

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INSTRUCTIONS ON ENVIRONMENTAL PROTECTION

Do not dispose of this product in the usual household garbage at the end of its life cycle; hand it over at a collection point for the recycling of electrical and electronic appliances. The symbol on the product, the instructions for use or the packing will inform about the methods for disposal.


The materials are recyclable as mentioned in its marking. By recycling, material recycling or other forms or re-utilization of old appliances, you are making an important contribution to protect our environment. Please inquire at the community administration for the authorized disposal location.

- ♦ Adjust the centering screws until the objective is perfectly centered.

IMPORTANT: before changing the lamp bulb or the fuse, make sure that the microscope is disconnected from the net. Check you are using the proper lamps, as other types can provoke malfunction.


Changing the lamp

- ♦ If the lamp is blown up, replace for another one, making sure it is a 12 V, 30 W Zuzi original one. Please take care and do not touch the lamp with bare hands. Do not use lamps with a higher power, as they could provoke over heating or any malfunction.
- ♦ For replacing the lamp, remove the precondenser (*Figure 5*)
- ♦ Replace the lamp, holding it firmly.

Changing the fuse

- ♦ If it is necessary to change the fuse loosen the protective lid that is at the back of the microscope over the socket; change the fuse and place it in the same position, putting the protective lid again. (*Figure 6*)

Cleaning

NOTE: Zuzi recommends using the microscope cleaner kit (code 89000001).

- ♦ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the microscope and produce an early ageing of the equipment.
- ♦ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

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
- ◆ Put the microscope in a horizontal plane stable table, having a safety area of at least 30 cm per side.
- ◆ Do not place the microscope near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.
- ◆ During operation dangerous materials such as flammable or pathological substances must be out of the working area.
- ◆ When you are not using the microscope for a long period of time please make sure it is unplugged from the net in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the microscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.


- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The microscope is supplied with a standard wire and it should be connected to a current wire provided with an earth wire, it should be handy to be disconnected in case of emergency.
- ◆ Never try to repair the microscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the microscope (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Zuzi Technical Assistance Department.

Centering the objectives

- ◆ Each objective can be centered individually in the nosepiece to align them perfectly with the light path and thus, avoid images distortions.
- ◆ For this purposes, insert the supplied keys in the centering screws (1.9) placed on each side of the objective.

Figure 4

- ◆ Adjust the condenser height (4.1) depending on the objective in use. The higher magnification of the objective the closer to the sample must be the condenser.
- ◆ Adjust the aperture of the field diaphragm (4.7) until it completely circumscribes the field of view.
- ◆ Adjust the aperture diaphragm (4.4) to the 70-80% of the numerical aperture of the objective in use.

Orthogonal observation

- ◆ The polarizer is located in the optical path so observations are performed under polarized light.
- ◆ Turn the polarizer (1.13) and the analyzer (1.6) until both are in position 0. Thus, both filters are in an orthogonal position (polarizer in E-W direction and analyzer in N-S direction).
- ◆ When using the 10x or a lower magnification objective remove the condensing lens (4.3) from the optical path. When using the 20x or a higher magnification objective put the condensing lens into the light path and raise the condenser to the uppermost position.
- ◆ Insert the complementary plates of gypsum, mica or quartz if necessary.

Conoscopic observation

- ◆ Conoscopic observation is performed with the 20x or a higher magnification objective.
- ◆ Put the condensing lens into the light path.
- ◆ Put the Bertrand lens (1.4) into the light path (position "B") and center it by using the centering screws (1.5).

Security

- ◆ The microscope must be used by previously qualified staff that knows how to operate it thanks to the user manual.

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
Figura 2


Figure 3

- ◆ Pre-focusing or focus limit control (*figure 3*): After focusing the sample with the low magnification objective, rotate the focus limit control (3.2) toward the front of the microscope to set the upper limit on the coarse adjustment movement. After changing samples or objectives, focusing is easily accomplished by rotating the coarse adjustment knob to reach the pre-focused position, then making fine adjustments with the fine adjustment knob. The fine adjustment is not affected by using the upper focus limit control.

Diopter compensation and interpupillary distance adjustment

- ◆ Diopter compensation: adjust to zero the dioptre compensation ring placed at the left eyepiece holder; then, looking through the right eyepiece, focus the sample specimen using the coarse and fine adjustment knobs. Once the sample is focused, look through the left eyepiece and adjust the dioptre compensation ring until obtaining a clear, sharp image.
- ◆ Interpupillary distance: Adjust the interpupillary distance by separating or joining the eyepieces until there is a total fusion of the two images.

Note: once the values of interpupillary distance and dioptre compensation are found it will be very useful to learn them by heart, mostly if the microscope is shared with other users. This is highly advisable to avoid repeating the adjustment each time the microscope is used.

Abbe condenser adjustment (Figure 4)

- 4.1 Condenser height adjustment knob
- 4.2 Condenser
- 4.3 Condensing lens lever
- 4.4 Aperture diaphragm

- 4.5 Condenser holding screw
- 4.6 Condenser centering screw
- 4.7 Pre-condenser with field diaphragm

◆ Put the binocular head and tighten firmly the holding screw (2.2) by using the supplied Allen key (Fig 2). Insert the eyepieces.

◆ Lower the stage by using the coarse adjustment knob (1.15). Screw the objectives on the nosepiece so as when it is turned clockwise the objectives follow an ascendant way (4x-10x-40x-60x).

◆ If necessary insert one of the compensator plates in the corresponding slot (1.7) in the intermediate tube. The plates are inserted at 45° from the preferential directions of vibration of the polarizer and analyzer.

◆ If wanted, put the attachable mechanical stage and fix it firmly.

◆ The microscope is supplied with a Schuko standard wire. Please insert the wire that feeds the AC electric current in the base of current 220 V, 50 Hz $\pm 10\%$ provided with earth wire and the other end to the microscope connector. The tension should be 220 V, 50 Hz $\pm 10\%$.

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220 V, 50 Hz $\pm 10\%$.

If you are not using the microscope for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

Operating mode

◆ Press the power on/off (1.17) and turn the adjusting knob (1.18) until the light intensity is adequate.

◆ Put the sample on the stage.

Focusing

◆ Using the objective with lesser magnification (4x), raise the stage by using the coarse adjustment knob (1.15) until the image of the sample appears. Then, move the fine adjustment knob (1.16) until the image is focused. By turning the nosepiece, put an objective with higher magnification (10x); it will only be necessary to adjust focusing with the fine adjustment knob.

◆ Focus tension adjustment (figure 3): Turning the friction control ring (3.1) toward the rear of the microscope increases the tension, and toward the front of the microscope loosens it. Avoid loosening it excessively since this could produce the fall of the stage together with problems on focusing accuracy.

Figure 3

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- Attachable mechanical stage
- Blue filter
- Lamp
- Fuse
- Objective centering keys
- Allen key
- Net wire
- Plastic cover
- User's manual
- Warranty certificate

We will only accept the equipment return within 15 days after delivery provided it comes in its original wrapping.

Installation


Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT.

- ♦ Please put the microscope on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the microscope near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc. Avoid the presence of flammable or toxic products in the working area.
- ♦ Put the intermediate tube and tighten firmly the holding screw (2.1) by using the supplied Allen key (*Fig. 2*).

2.1 Intermediate tube holding screw

2.2 Head holding screw

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Code	50171146
Model	146P
Optics	Infinity corrected
Head	Binocular Siedentopf-type; 30° inclined and 360° rotary Interpupillary distance 55-75 mm; Diopter compensation
Eyepieces	2 eyepieces WF10x/22 mm 1 eyepieces WF10x/20 mm graduated
Nosepiece	Quadruple: 4x, ∞ (A.N.: 0.10); 10x, ∞ (A.N.: 0.25); 40x (R), ∞ (A.N.: 0.65); 60x (R), ∞ (A.N.: 0.85)
Objectives	Plan achromatic; can be centered individually
Polarizer	Rotary 360°
Analyzer	Rotary 90°
Bertrand lens	Removable from optical path
Compensator plates	Gypsum (1 λ), mica (1/4 λ) and quartz wedge.
Stage	Round (160 mm Ø), rotary 360° and graduated (resolution 0.1°) Fixed clamps or mechanical stage with movement range 35x25 mm
Condenser	Abbe (N.A.: 1.25) with condensing lens removable from the optical path
Diaphragms	Field and aperture
Lamp	Halogen 12 V / 30 W. CA 85 V – 230 V.

4. INSTALLATION / SETTING UP

Unwrap the microscope, take off the involving plastic and take off the polyspan protection in which it comes fitted. Take off all the protective items and, without connecting the microscope to the net, make sure that it does not present any damage because of the shipment. In case the microscope presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- Binocular head
- Intermediate tube
- 2 Eyepieces WF10x/22 mm and 1 eyepiece WF10x/20 mm graduated
- Plan achromatic objectives: 4x, 10x, 40x and 60x
- Compensator plates: gypsum (1 λ), mica (1/4 λ) and quartz wedge

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
2. DESCRIPTION

- | | |
|---|------------------------------|
| 1.1 Binocular head | 1.10 Rotary stage |
| 1.2 Eyepieces | 1.11 Stage fixing screw |
| 1.3 Intermediate tube | 1.12 Condenser |
| 1.4 Bertrand lens | 1.13 Polarizer |
| 1.5 Bertrand lens centering screws | 1.14 Pre-condenser |
| 1.6 Analyzer | 1.15 Coarse adjustment knob |
| 1.7 Slot for inserting compensator plates | 1.16 Fine adjustment knob |
| 1.8 Nosepiece | 1.17 Power on/off |
| 1.9 objective centering screws | 1.18 Light intensity control |

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INDEX OF CONTENTS

1. USES OF THE INSTRUMENT	14
2. DESCRIPTION	15
3. TECHNICAL SPECIFICATIONS	16
4. INSTALLATION / SETTING UP	16
5. MAINTENANCE AND CLEANING	21
ANNEX I: CE CERTIFICATE	24

1. USES OF THE INSTRUMENT

Zuzi microscope 146P equipped with infinity corrected optics offers excellent optical features for observing samples under polarizing light, being of application in petrography, mineralogy and materials chemistry as well as in histology and anatomical pathology for studying birefringent structures (crystalline structures, pigments, lipids, amyloid deposits, etc.).

For this purpose, the microscope is provided with a polarizing filter under the stage that only transmits light in N-S direction, and an analyzing filter, located over the stage, that transmits light in E-W direction. When both filters are in orthogonal position (at an angle of 90°), the light generated by the polarizer is extinguished and nothing is observed; however, if the sample is birefringent, the structure will light up on a dark field.

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Zuzi polarizing microscope 146P. We highly recommend looking after this equipment according to what is stated in this manual.

Zuzi develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ♦ This manual is inseparable from the Zuzi polarizing microscope 146P, so it should be available for all the users of this equipment.
- ♦ You should carefully handle the microscope avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ♦ Never dismantle the different pieces of the microscope to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the microscope.
- ♦ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ♦ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Zuzi Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0.30 Euros/min).
- ♦ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ♦ Overhaul is not covered by the microscope warranty.
- ♦ Operations made by non-qualified staff will automatically produce a loss of the microscope warranty.
- ♦ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ♦ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Zuzi Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ♦ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ♦ Manufacturer reserves the right to modify or improve the manual or equipment.